

EST.
2015




SHOE
MFG

ANODYNE


PREMIUM QUALITY

Anodyne Scanner App Instructions



Once you've downloaded the Anodyne App and attached your Structure camera, you're ready to go! Below are step by step instructions on how to properly submit an order via the Anodyne App.

1. Click on the Anodyne App Icon to launch the scanner application.




2. If this is the first time opening the App, you will be prompted to enter new login credentials, these should have been sent to you from the Anodyne team already, if not, please reach out to customer service, 844-637-4637, and credentials will be provided.

A screenshot of the Anodyne app login screen. The screen is white with a grey header bar containing a close button (X) and the word "LOGIN" in orange. Below the header, the text "HAVE AN ACCOUNT?" is displayed in orange, followed by the instruction "Log in to your account by entering your credentials below" in grey. There are two input fields: "Email" and "Password". Below the password field, there is a link that says "Forgot password? Or need to create an account? Just give us a call at 1-844-637-4637 and we'll be happy to help!". At the bottom of the screen, there is a large orange button labeled "LOG IN".

3. Select the “User Info” button located on the bottom right of the screen.



4. Now that you’re in the “User Info” menu, you will be able to see your pre-loaded account information, which is tied directly to your Customer Account with Anodyne. If this information needs to be changed, please reach out to our Customer Support team, 844-637-4637



5. You're now ready to scan a patient. Select the "New Patient" button located on the bottom left of the screen, and then click the type of inserts you would like to order.




6. After you select which style of inserts you like, you will be taken to the "Order Form." At the very minimum, fill out all of the fields that are highlighted in red. Those fields are required to be filled in before you can move on to the next step. After this has been completed, click "Next," which will take you to the "Accommodations" section.

The screenshot shows the "ORDER FORM" page. At the top right is a "SAVE" button with a download icon. Below the header is the "ORDER DETAIL" section, which contains four red-outlined input fields: "Patient Name", "Gender" (dropdown), "Activity" (dropdown), and "Weight" (dropdown). Below these is a "PO Number" input field. The "INSERT DETAIL" section follows, with two rows of dropdown menus:

- Row 1: Insert Type-No. 3 Custom Accommodative Inserts, Quantity-3, Pairs
- Row 2: Base-45 EVA, Top Cover-Diabetic, Length-Full, Midlayer

A red arrow points from the bottom left towards the "NEXT >" button, which is highlighted in orange. The navigation bar at the bottom is identical to the one in the previous screenshot, containing: NEW PATIENT, SAVED PATIENTS, SUPPORT, and USER INFO.

7. On the “Accommodations” section, you can add any special accommodations or instructions for the inserts, which best fit your patient’s needs. After you’ve added anything, hit Next.


ORDER FORM SAVE 





ACCOMMODATIONS

Arch Height	High Lateral Flange	Offload/Sweetspot Metatarsal
Heel Cup	Medial Arch Flange	Left: <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5
Met Pad	Offload Heel	Right: <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5
Met Bar	Offload Arch	Offload/Sweetspot Distal
Heel Raise	Rear Foot Posting	Left: <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5
Left	Left	Right: <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5
Right	Right	Charcot Deformity
		Additional Comments
		<div></div>


< BACK

NEXT >




 NEW PATIENT  SAVED PATIENTS  SUPPORT  USER INFO

8. You will then be taken to the Order Confirmation section. Here you will see a summary of all the patient information entered so far. Please make sure to double check your order details and select an address for the order to ship to. This is done, by clicking the pull down menu on the top, and selected the desired address from the options available.

ORDER OVERVIEW SAVE 

ORDER CONFIRMATION



MILWAUKEE-5050 S 2ND STREET

Patient Name: Michael P

PO:

Shoes: M038-40-M-075


Inserts: No. 3 Custom Accommodative Inserts


Insert Qty: 3


Account Number: 5050

Account Name: michaelp@anodyneshoes.com

Email: michaelp@anodyneshoes.com





 SCAN LEFT FOOT

 SCAN RIGHT FOOT

 ADD PHOTOS


< BACK

UPLOAD

 NEW PATIENT  SAVED PATIENTS  SUPPORT  USER INFO

9. If the desired ship to location is not listed, select the “Add Custom Address” option, and enter the address. This will save to your addresses for future use as well

ORDER OVERVIEW


SAVE 

ORDER CONFIRMATION


ADD CUSTOM ADDRESS

Patient Name: Michael P
PO:
Shoes: M038-40-M-075
Inserts: No. 3 Custom Accommodative Inserts
Insert Qty: 3


Account Number: 5050
Account Name: michaelp@anodyneshoes.com
Email: michaelp@anodyneshoes.com



SCAN LEFT FOOT




SCAN RIGHT FOOT





ADD PHOTOS


< BACK

UPLOAD


 NEW PATIENT

 SAVED PATIENTS

 SUPPORT

 USER INFO

ORDER OVERVIEW

SAVE 

ORDER CONFIRMATION

ADD CUS

Name

Address Label

Address 1

Address 2

City

State

Zip Code

CANCEL

CREATE ADDRESS


SCAN LEFT FOOT


SCAN RIGHT FOOT


ADD PHOTOS


< BACK

UPLOAD

 NEW PATIENT


 SAVED PATIENTS

 SUPPORT

 USER INFO

10. Now you are ready to start scanning. Select the “Scan Left” or “Scan Right” icon.

ORDER OVERVIEW

SAVE 

ORDER CONFIRMATION

MILWAUKEE-5050 S 2ND STREET

ADD CUSTOM ADDRESS

SET TO DEFAULT

Address Name: Anodyne

Address 1: 5050 S 2ND STREET

Address 2: - None -

City: MILWAUKEE

State/Province: WI

Zip Code: 53207

Patient Name: Michael P

PO:

Shoes: M038-40-M-075


Inserts: No. 3 Custom Accommodative Inserts

Insert Qty: 3


Account Number: 5050

Account Name: michaelp@anodyneshoes.com


Email: michaelp@anodyneshoes.com



SCAN LEFT FOOT




SCAN RIGHT FOOT





ADD PHOTOS


< BACK

UPLOAD


 NEW PATIENT

 SAVED PATIENTS


 SUPPORT

 USER INFO


11. With the patient properly positioned, you’re ready to begin your scan. Once the scanner is mapping the blue/green to the bottom of the foot, click the “Scan” button. Zooming in or out by pinching or spreading your fingers will allow you to change the size of the 3D box. You’ll have to do this, in addition to moving closer or further from the patient in order to optimize the mapping.



12. After you select “Scan,” the 3D image will begin rendering. Slowly move the camera around the foot to make sure that all of the anatomy is scanned. Once you’re satisfied with the scan, click the “Done” icon.




13. In the “Scan Result Viewer” you can inspect the scan by moving the 3D image around. Select “Back”, located in the top left corner, if you are unhappy with the image. Otherwise, click the “Approve” icon, in the top right corner, to proceed.




14. You'll then be taken back to the "Patient Form" where you will scan the other foot and repeat steps 9 through 11.

ORDER OVERVIEW

SAVE 

ORDER CONFIRMATION

MILWAUKEE-5050 S 2ND STREET 


ADD CUSTOM ADDRESS

SET TO DEFAULT


Patient Name: Michael P
PO:
Shoes: M038-40-M-075
Inserts: No. 3 Custom Accommodative Inserts
Insert Qty: 3

Address Name: Anodyne
Address 1: 5050 S 2ND STREET
Address 2: - None -
City: MILWAUKEE
State/Province: WI
Zip Code: 53207


Account Number: 5050
Account Name: michaelp@anodyneshoes.com
Email: michaelp@anodyneshoes.com



RE-SCAN
LEFT FOOT




SCAN RIGHT FOOT





ADD PHOTOS


< BACK

UPLOAD

 NEW PATIENT


 SAVED PATIENTS

 SUPPORT


 USER INFO

15. In addition to taking the scans, our App also gives you the option of sending us photos of the patient's feet. To add Photos, click the "Add Photos" icon.

ORDER OVERVIEW

SAVE 

ORDER CONFIRMATION

MILWAUKEE-5050 S 2ND STREET 


ADD CUSTOM ADDRESS

SET TO DEFAULT


Patient Name: Michael P
PO:
Shoes: M038-40-M-075
Inserts: No. 3 Custom Accommodative Inserts
Insert Qty: 3

Address Name: Anodyne
Address 1: 5050 S 2ND STREET
Address 2: - None -
City: MILWAUKEE
State/Province: WI
Zip Code: 53207


Account Number: 5050
Account Name: michaelp@anodyneshoes.com
Email: michaelp@anodyneshoes.com



RE-SCAN
LEFT FOOT




RE-SCAN
RIGHT FOOT





ADD PHOTOS


< BACK

UPLOAD


 NEW PATIENT

 SAVED PATIENTS


 SUPPORT

 USER INFO

16. Next, select the Camera Icon to take a new photo




17. Or you can choose to use a previously saved picture. Click "SET", and you will have the ability to select a photo that is saved to your tablet



18. Once you've selected or taken your photo you will have the ability to use your finger or stylus to draw on it by clicking "Edit". This helps to highlight any sort of specific complications, wounds, and/or call-outs. Click "Save and Exit" when you are satisfied.




19. Once you're satisfied with the picture(s) that you've added, select "Back" in the top left corner of the screen, and you will be taken back to the "Patient Form."



20. After both feet have been successfully scanned, the order form has been filled out, and/or pictures have been added, you're ready to send everything over to us. Review all of the information one last time and click the "Upload" icon to send the files.

ORDER OVERVIEW

SAVE 

ORDER CONFIRMATION

MILWAUKEE-5050 S 2ND STREET

ADD CUSTOM ADDRESS

SET TO DEFAULT

Address Name: Anodyne

Address 1: 5050 S 2ND STREET

Address 2: - None -

City: MILWAUKEE

State/Province: WI

Zip Code: 53207

Patient Name: Michael P

PO:

Shoes: M038-40-M-075


Inserts: No. 3 Custom Accommodative Inserts

Insert Qty: 3


Account Number: 5050

Account Name: michaelp@anodyneshoes.com


Email: michaelp@anodyneshoes.com



RE-SCAN LEFT FOOT




RE-SCAN RIGHT FOOT





ADD PHOTOS


< BACK

UPLOAD


 NEW PATIENT

 SAVED PATIENTS


 SUPPORT

 USER INFO


21. After you click the "Upload" icon, you'll receive a pop-up notifying you that your order was sent successfully.




COMFORT.
ALL DAY. EVERY DAY.




No. 3
Custom Accommodative Inserts
A5514 Reviewed



No. 4
Partial Foot Toe Filler
L5000 Reviewed



No. 5
Semi-Rigid Orthotics





No. 6
Custom Accommodative Cork Inserts
A5514 Reviewed


SUBMIT SUCCESSFUL


Order created

OK

 NEW PATIENT

 SAVED PATIENTS

 SUPPORT

 USER INFO

22. If you would like to revisit or check scans/orders for previous patients, select the “Saved Patients” option on the home screen.



23. Here you can see when previous orders were sent.

SAVED PATIENTS				
SORT BY DATE		SORT BY NAME		
Not Uploaded	Billy	View Patient	DELETE	
Not Uploaded	Test	View Patient	DELETE	
07/16/2020	Michael P	View Patient	DELETE	
07/09/2020	Bobby	View Patient	DELETE	
06/09/2020	Test	View Patient	DELETE	

NEW PATIENT | **SAVED PATIENTS** | SUPPORT | USER INFO